

# SESLHD PROCEDURE COVER SHEET



**Health**  
South Eastern Sydney  
Local Health District

<b>NAME OF DOCUMENT</b>	Sterilisation: Routine Monitoring of Cleaning, Disinfection and Sterilisation Processes
<b>TYPE OF DOCUMENT</b>	Procedure
<b>DOCUMENT NUMBER</b>	SESLHDPR/548
<b>DATE OF PUBLICATION</b>	May 2020
<b>RISK RATING</b>	Low
<b>LEVEL OF EVIDENCE</b>	National Safety Quality Health Service Standards: Standard 3 – Preventing and Controlling Healthcare Associated Infection  Australian/New Zealand Standard AS/NZS 4187:2014
<b>REVIEW DATE</b>	May 2025
<b>FORMER REFERENCE(S)</b>	N/A
<b>EXECUTIVE SPONSOR or EXECUTIVE CLINICAL SPONSOR</b>	Director of Clinical Governance and Medical Services
<b>AUTHOR</b>	SESLHD Sterilising Services Working Party (SSWG)
<b>POSITION RESPONSIBLE FOR THE DOCUMENT</b>	Manager Sterilising Services, The Sutherland Hospital Karolina.Tipevska@health.nsw.gov.au
<b>KEY TERMS</b>	Reusable Medical Device (RMD), Routine Monitoring, Validation, High Level Disinfection (HLD), Preformed Sterile Barrier System (PSBS), Process Challenge Device (PCD), Sterile Barrier System (SBS)
<b>SUMMARY</b>	To ensure processes undertaken by the Sterilising Services Departments consistently produce a result that meets predetermined specifications.

## **COMPLIANCE WITH THIS DOCUMENT IS MANDATORY**

**This Procedure is intellectual property of South Eastern Sydney Local Health District.  
Procedure content cannot be duplicated.**

Feedback about this document can be sent to [seslhd-executiveservices@health.nsw.gov.au](mailto:seslhd-executiveservices@health.nsw.gov.au)

# SESLHD PROCEDURE

## Sterilisation: Routine Monitoring of Cleaning, Disinfection and Sterilisation Processes

SESLHDPR/548

### 1. POLICY STATEMENT

The purpose of routine monitoring and control is to provide evidence that the specified and validated cleaning, disinfection, packaging and sterilisation processes for an RMD have been achieved.

### 2. BACKGROUND

Installation Qualification (IQ), Operational Qualification (OQ), Performance Qualification (PQ) and Validation has been performed on Washer/Disinfectors, Sterilisers and the associated equipment.

### 3. RESPONSIBILITIES

#### 3.1 Employees will:

- Comply with the requirements of this procedure
- Report any non-compliance to the Sterilising Services Manager.

#### 3.2 Line Managers will:

- Support compliance

#### 3.3 District Managers/ Service Managers will:

- Implement the requirements of this procedure to assure the quality and safety of reprocessed RMDs
- Ensure routine monitoring is performed and recorder.

#### 3.4 Medical staff will:

- NIL

### 4. PROCEDURE

#### 4.1 Routine Monitoring and Control of the Cleaning Process:

##### 4.1.2 Manual Cleaning

The outcome of the manual cleaning shall be checked at completion of process by visual inspection.

##### 4.1.3 Washer/Disinfectors Employing Thermal Disinfection

Washer/Disinfectors shall be tested each day, at the beginning of the day, to ensure that is functioning as intended.

The following parameters shall be verified at completion of each cycle:

- Visual inspection
- Result of process challenge devices
- Process/cycle parameters such as: temperature and exposure time.

## SESLHD PROCEDURE

### Sterilisation: Routine Monitoring of Cleaning, Disinfection and Sterilisation Processes

SESLHDPR/548

#### 4.1.4 Ultrasonic/Irrigator

Following testing shall be performed to ensure Ultrasonic/Irrigator is functioning as intended:

- Cleaning efficacy as per manufacturer's instructions
- Performance of individual transducers as per manufacturer's instructions

#### 4.1.5 Drying Cabinets:

Temperature shall be checked and recorded on daily basis, ensuring within parameters.

### 4.2 Routine Monitoring and Control of Processes Employing Chemical Disinfection for Thermolabile RMDs

4.2.1 Routine monitoring of the HLD Process as per manufacturer's instructions, including but not limited to:

- Performance of disinfection equipment-water pressure, flow and temperature
- Chemical agent concentration
- Contact time
- Rinse time/volume.

### 4.3 Microbiological Surveillance of Flexible Endoscopes with channels

- Flexible Endoscopes shall undergo microbiological surveillance in accordance with HSO Procedures and GENCA Guidelines.

### 4.4 Routine Monitoring and Control of Packaging Processes

- SBS
- PSBS
- Reusable rigid sterilising containers
- Heat sealer temperature shall be recorded daily and visual checks shall be made immediately after each episode of sealing
- One or more samples shall be checked daily to establish seal integrity before and after exposure to a sterilisation process.

### 4.5 Routine Monitoring and Control of Sterilising Processes

Routine monitoring and control shall demonstrate that the specified and validated sterilisation processes for an RMD have been delivered to that RMD

#### 4.5.1

- Sterilising equipment shall be checked according to Table 8.2 to ensure that it is functioning as intended each day
- Process records shall be checked at the completion of each cycle to verify that the process was delivered in accordance with specifications
- Biological Indicators (BI), Chemical Indicators (CI) and Process Challenge Devices (PCD) may be used as an additional method demonstrating sterilisation processes were delivered in accordance with the validated specifications.

# SESLHD PROCEDURE

## Sterilisation: Routine Monitoring of Cleaning, Disinfection and Sterilisation Processes

**SESLHDPR/548**

### 4.5.2 Low temperature Sterilisation Systems

- Performance tests shall be conducted in accordance with the manufacturer’s instructions for use.

### 4.5.3 Moist Heat

- Daily air removal and steam penetration test (Bowie and Dick Type Test) shall be performed on steam sterilisers
- Leak rate/vacuum test and air detector tests shall be performed weekly on steam sterilisers with an air detector, daily for steam sterilisers without air detector
- Requirements for routine monitoring and control of sterilising equipment as per Table 8.2.

## 5. DOCUMENTATION

- Specific HSO Procedures

## 6. AUDIT

- Reusable Medical Devices (RMDs) - Packing and Sterilisation of - Daily Audit
- Reusable Medical Devices (RMDs) - Mechanical cleaning of - Daily Audit
- Reusable Medical Devices (RMDs) - Manual cleaning of - Daily Audit

## 7. REFERENCES

- AS/NZS 4187:2014 Reprocessing of reusable medical devices in health services organizations
- ISO 11607 Packing for terminally sterilised medical devices
- ISO11135 Sterilisation of health care products-Biological Indicators (series)
- ISO11140 Sterilisation of health care products-Chemical Indicators (series)
- ISO15882 Sterilisation of health care products-Chemical indicators-Guidance for selection, use and interpretation of results
- ISO17665 Sterilisation of health care products-Moist heat
- ISO17665-1 Part 1: Requirements for the development , validation and routine control of a sterilisation process for medical devices
- ISO25424 Sterilisation of health care products-Low temperature steam and formaldehyde - Requirements for development, validation and routine control of a sterilisation process for medical devices
- ISO14180 Sterilisers for Medical Purposes-Low Temperature Steam Formaldehyde Sterilisers – Requirements and Testing
- EN 285 Sterilisation-Steam Sterilisers - Large Sterilisers

## 8. REVISION AND APPROVAL HISTORY

Date	Revision No.	Author and Approval
November 2016		SESLHD Sterilising Resource Group
November 2016	1	Approved by Executive Sponsor for Draft for Comment

**SESLHD PROCEDURE****Sterilisation: Routine Monitoring of Cleaning,  
Disinfection and Sterilisation Processes****SESLHDPR/548**

December 2016	1	Approved by Clinical and Quality Council
February 2020	2	SESLHD Sterilising Working Party (SSWG) conducted a minor review. Reworded to more accurately reflect the current standards. References updated and listed.
May 2020	2	Approved by Executive Sponsor. Published by Executive Services.